The new Dimension for the Laboratory

Coagulation Analyzer

- **Medium laboratories**
- **Doctors office**
- Intensive care
- Research facility

REF: 80 700 00



- with sunlight reduction
- Clotting, chromogenic or immunturbidimetric assays
- **Quantitative D-Dimer**
- Micro volumes (total 75 μL)
- Small dimension and weight
- Easy operation
- 6 languages on board
- Highly reliable & long-living
- Autosense optics to eliminate interferences like Bilirubin, Hemoglobin, sunlight
- Approved clotting algorithm with biphasic waveform analysis
- Automatic start at reagent addition
- 5 point calibration curve for each test
- 2 Stop-watch function onboard
- **Double determination**
- Magnetic stirrer for reagent
- **Profile testing like PT+APTT**
- Optional external thermal printer
- Optional barcode scanner
- Optional automatic start pipette
- RS232 interface for data research or firmware upgrade
- **Linking to Laboratory Information** Systems such as TECO TECAM















Available Tests

РΤ Prothrombin Time

Activated Partial Thromboplastin Time APTT

FIB Fibrinogen (Clauss and derived)

Thrombin Time TT VT Venom time

II,V,VII,VIII,IX,X,XI,XIIFAC

Antithrombin AT3 **HEP** Heparin PC Protein C PS Protein S PLG Plasminogen A2AP a2-Antiplasmin

PCA Protein-C resistance (automatic ratio) Screen & confirm (automatic ratio) LA

ECAT Ecarin Chromogenic Assay For Thrombin ECAH Ecarin Chromogenic Assay For Hirudin

quantitative D-Dimer DD

Available dimensions

Seconds, %, Ratio, INR, mg/L, mg/dL, g/L, U/mL, ng/mL, µg/mL

Development, Manufacturing and Sales:

TECO Medical Instruments, Production + Trading GmbH

Dieselstrasse 1 D-84088 Neufahrn NB

Internet: http://www.teco-gmbh.com Email: <u>info@teco-gmbh.c</u> Phone: +49 8773 910010 info@teco-gmbh.com Fax: +49 8773 910011





Distributor:

Technical Data

Dimensions:

LxWxH: 230x200x110mm

Weight: 1,40kg

Power Supply (external):

Input: 100-240V~, 47-63Hz

Power: 42W

Approvals: CE, EN60601-1, UL

Incubation Block:

Temperature: $37,0^{\circ}\text{C} \pm 0,2^{\circ}\text{C}$ Samples: 12 cuvettes

Reagent: 1 vial Ø30,0mm, stirred

1 vial Ø22,5mm 1 vial Ø16mm

Different adapters available

Measurement:

Optic: 2 channel with AutoSense
Methods: Turbidity, Kinetic, Endpoint
Algorithm: based on optical density
Processing: 18Bit ADC + Digital Filter
Autosense: LED intensity controlling
Starting: manual, optical, autopipette

Performance:

Precision: CV < 5% Throughput PT ~70/h APTT ~30/h

FIB ~50/h

Optional Accessory

• Autopipette:

Volume: 4-stage 25/50/100/200µL

Test-Start: electronical

Barcode Scanner:

Interface: serial, RS232
Type: CCD Handheld
Power: powered over RS232

Printer:

Interface: serial, RS232

Type: external thermal printer

TECAM SMART for MS Windows XP, Vista

(Laboratory Information System)

A perfect and user friendly tool to combine laboratory data management, quality control and research purpose in one.

Intuitive handling

One click principle for data collection, filtering, reporting or QC

Interface

Serial RS232 to analyzer. The results will be transferred including the optical reaction curve

Database

Manage over 1Mio results and patient information in one database

Stored Information

Date, Patient (Id, name, sex, birth), Test, Lot, Results, Flags, Error, System, Range, Comment and Reaction curve

Complete Traceability by one click

Review for Date, PID, Test, QC by one click

 Statistical analyse for QC or patients Westgard Multirules Analysis Levey Jenning diagram

• Flexible Multifilter

Combine every data field. Predefined timefilters (eg. today, last 30 days, morning, afternoon etc.)

Reporting

make reports by one click. Use multifilter or presets to generate specific reports (eg. show results of specific patient of all PT within the last 30 days)

Research

Export data, grids, charts, curves to PDF, XLS, TXT or HTML file

Security

Secured access to critical menus

- Online update
- Multi-Lingual

Standard package

1 Pc Coatron M2 1 Pc Power Supply 25 Pcs Double cuvettes 1 Pc Mains cable 5 Pcs Reagent container 1 Pc Magnetic stirrer 1 Pc Printer cable Reagent adaptor 1 Pc 1 Pc Pipette fastener 1 Pc Dust cover 1 Pc **Operation Manual** 1 Pc Warranty card 1 Pc Installation report 1 Pc Service report